



SurgiCare Recessed Sealed 2x4

2x4 Sealed Supplemental Surgical Luminaire, SSTR24 Series

2x4 Sealed Supplemental Surgical Luminaire, SSTR24 Series

The SSTR Series is designed for the illumination, electrical shielding and infection control requirements of the surgical suite. The asymmetric/symmetric light distribution provides the ambient illumination while contributing supplemental illumination for the surgical field. SSTR luminaires are available in a variety of lamping options. The housing is fully sealed to minimize potential room contamination and simplify cleaning. Meets MIL-STD 461F, wet location listed and is IP65 rated.

Product data

General Information		Shape	
Lamp Envelope Code	T-8	Rectangular	
Lamp Source Code	Fluorescent		
Lamp Watts	32		
Number of Lamps	6		
Photometric Distribution Type	Direct		
Operating and Electrical		Product Data	
Voltages	120/277 V	Order product name	2x4 Sealed Supplemental Surgical Luminaire, SSTR24 Series
		Order code	SSTR(F/G)24-6-T8-ASM
		Material Nr. (12NC)	ALK00004173
		Catalog Number	SSTR(F/G)24-6-T8-ASM
		Catalog Number Description	2x4 Sealed Supplemental Surgical Luminaire, SSTR24 Series
Mechanical and Housing			
Fixture Size	2 ft x 4 ft		