



# Product Description

## Medical Therapy Jaundice TL/TL-D

Tubular or compact single-ended low-pressure mercury vapour fluorescent lamps

### Benefits

- Optimal spectrum for photo-oxidative process to convert unconjugated bilirubin into a watersoluble form

### Features

- Emit radiation for therapeutic and other applications

### Application


- Medical treatment of jaundice in new-born babies (hyperbilirubinaemia)

### Versions




# Medical Therapy Jaundice TL/TL-D

## Dimensional drawing



Product	D (max)	A (max)	B (max)	B (min)	C (max)
TL-D 18W/52	28 mm	589.8 mm	596.9 mm	594.5 mm	604 mm



Product	D (max)	A (max)	B (max)	B (min)	C (max)
TL 20W/52 SLV/25	40.5 mm	589.8 mm	596.9 mm	594.5 mm	604 mm

