



# Bio-Seal - Superior and practical

## BioSeal Recessed Cleanroom

Bio-Seal is designed for the most demanding biological, chemical, immunological, or food chemistry laboratories. It is suitable for hot, cold, and sealed locations with listings for cleanroom class 10-10,000, NSF, and IP66.

### Benefits

- Meets multiple needs
- Ease of maintenance

### Features

- Class 10 to 10,000 cleanrooms
- IP66
- Suitable for wet locations
- Available in a variety of sizes and materials

### Application

- Cleanrooms
- Food Processing
- Hazardous Location

## BioSeal Recessed Cleanroom

### Versions

