



Meets all illumination needs for a surgical suite

SurgiCare Recessed 2x4

The STR Series is designed for the illumination, electrical shielding and infection control requirements of the surgical suite. The asymmetric light distribution provides the ambient illumination while contributing supplemental illumination for the surgical field. STR luminaires are available in a variety of lamping options. Meets MIL-STD 461F.

Benefits

- One piece door frame
- antimicrobial pain finish
- smooth, room side lens

Features

- Meets MIL-STD 461F

Application

- Hospitals and Care Centers

SurgiCare Recessed 2x4

Versions

