



# LITEGUARD

Minimizing the risk of infection in healthcare environments is a top priority. The Liteguard helps protect patients, staff, visitors from undesirable pathogens and allergens.



The Liteguard complements almost any recessed luminaire of 2'x2', 2'x4' and 1'x4'.



**PHILIPS**  
**LIGHTOLIER®**

# LITEGUARD LG

## Protection for Recessed Luminaires



Neoprene gasket seals the frame to the T-Bar, preventing the collection of unwanted particles on the louvers or inside the center basket of the luminaire.

### Features

- Lens for recessed luminaires adding protection against pathogens, dust and bugs and makes cleaning easy.
- Accessory with low upfront price and reduces cleaning time and costs
- Fits on almost all 2'x2', 2'x4' and 1'x4' fully recessed luminaires including troffers, parabolics and Alter-types.
- Light output reduction of only  $\pm 5\%$  with the clear lenses on Alter and parabolics luminaires.
- Post-painted, matte white aluminum frame of 1 1/2" thick blending perfectly with T-Bars and white fixtures.
- Complete easy-to-install kit including T-bar supports and screws. Contractor responsible for sealing fixture to ceiling.
- Can be used on drywall ceilings with flange kits (order separately according to fixture size FK91X4, FK92X2 or FK92X4)

### Ordering Guide

Example : LG2X2GI

**LG**

**Family**  
Liteguard

**2X2**

**Luminaire Size**  
**2X2** : 2'x2'  
**2X4** : 2'x4'  
**1X4** : 1'x4'

**GI**

**Lens**  
**GI** : Clear Acrylic  
**PX** : Clear Polycarbonate



**GI** : Clear Acrylic Lens

**PX** : Clear Polycarbonate Lens

Contact your Canlyte representative for more information.

1-800-668-2770 • [www.canlyte.com](http://www.canlyte.com)

We reserve the right to change details of design, materials and finish.

© 2009 Philips Group • CG307E