



Image  
not yet  
available

# Metronomis Torino

## ZDS500 LO

### METRONOMIS ACC

Metronomis is a family of ten highly distinctive luminaires and dedicated supports (masts and brackets), which can be used with a wide variety of optics to provide architects and lighting designers with the total solution needed to create a unified, consistent lighting design that still reflects the differences in urban culture and history. In its daytime appearance, Metronomis integrates harmoniously with any street scene, while its night-time appearance is both functional – guiding people and traffic safely through the city – and decorative, creating a pleasant ambience in city squares and pedestrian precincts. As the current range is based on conventional light sources, four highly innovative, LED-dedicated designs have been developed, echoing the original silhouette but retaining flexibility and modularity as key characteristics. The Metronomis LED range offers a number of new ways to personalize your urban projects. Please check out the Metronomis LED luminaire variants.

### Product data

| General Information |                         |                                 |               |
|---------------------|-------------------------|---------------------------------|---------------|
| Product family code | ZDS500 [METRONOMIS ACC] | Material Nr. (12NC)             | 910502250018  |
|                     |                         | Numerator - Quantity Per Pack   | 1             |
|                     |                         | EAN/UPC - Product/Case          | 8717943952816 |
|                     |                         | Numerator - Packs per outer box | 1             |
|                     |                         | EAN/UPC - Case                  | 8717943952816 |
| Product Data        |                         |                                 |               |
| Order product name  | ZDS500 LO               |                                 |               |
| Full product name   | ZDS500 LO               |                                 |               |
| Full product code   | 871794395281600         |                                 |               |
| Order code          | 910502250018            |                                 |               |