



# StyliD PerfectBeam

## ST741Z LBSP

StyliD PerfectBeam has been developed specifically for museum applications. It strikes the perfect balance between presenting art and preserving it, between ever-changing exhibitions and the need for low maintenance and energy efficiency. The often conflicting demands for good visibility and conservation are now being met with LED lighting that emits no harmful UV or IR radiation. StyliD PerfectBeam presents works of art in their true colors, subtly attracting the viewer's eye, making it easy to see, interpret and appreciate masterpieces down to the smallest detail. It enables you to use light to evoke emotions and create an unforgettable museum experience. Its modular design, variable beam, beam shaping accessories and onboard dimming ensure maximum flexibility to cope with changing exhibitions and lighting needs, now and in the future.

### Product data

Product Data			
Order product name	ST741Z LBSP	Material Nr. (12NC)	910504801189
Full product name	ST741Z LBSP	Numerator - Quantity Per Pack	1
Full product code	871829129783300	Numerator - Packs per outer box	1
Order code	910504801189		