



# Latina FBH024

## GBH024 RG

Latina is a recessed downlight for compact fluorescent PL-C lamps, with spring fasteners for easy and quick installation. Square and round models are available. Latina comes in ready-to-install kit versions, including lamps and pre-connected separate gearbox. The round models are available in glass versions (frosted glass or clear decorative glass) or in open versions. For the open versions, the same glass designs are available as accessories (to be ordered separately). The square versions are always with glass.

### Product data

Product Data	
Full product code	871794388532999
Order product name	GBH024 RG
Order code	910930025318
Numerator - Quantity Per Pack	1
Numerator - Packs per outer box	20

Material Nr. (12NC)	910930025318
Full product name	GBH024 RG
EAN/UPC - Case	8717943886760
EAN/UPC - Product/Case	8717943885329