



# ClassicStreet brackets

## JGB794 MBP-T 750 76 BK

CLASSICSTREET BRACKETS SUSPENDED, Pole-mounting bracket 2 arms 180°, 750 mm, Black

Recalling the iconic beauty of 20th-century street lanterns, ClassicStreet combines elegant design with the subtle details of heritage luminaires while meeting today's urban lighting requirements. Though echoing the past, the ClassicStreet luminaire is very much shaped for the future. Its gently curved design, together with dedicated brackets and poles, creates a contemporary look that beautifies our cities and creates distinctive and warm night- and daytime ambiances. Not only in historical and traditional architectural settings, but also in more modern urban environments too. Designed around its LED engine, ClassicStreet delivers outstanding levels of lighting and energy performance.

### Product data

General Information		Full product name	JGB794 MBP-T 750 76 BK
Bracket type	Pole-mounting bracket 2 arms 180°	Full product code	871794321652900
Bracket color	Black	Order code	21652900
Product family code	JGB794 [CLASSICSTREET BRACKETS SUSPENDED]	Material Nr. (12NC)	912300023375
Value ladder	Specification	Numerator - Quantity Per Pack	1
CE mark	Yes	EAN/UPC - Product/Case	8717943216529
Mechanical and Housing		Numerator - Packs per outer box	1
Bracket length	750 mm	EAN/UPC - Case	8717943216529
Product Data			
Order product name	JGB794 MBP-T 750 76 BK		

## ClassicStreet brackets

### Dimensional drawing

