



ClassicStreet brackets

JGB795 MBP-T 750 76 BK

CLASSICSTREET BRACKETS POST TOP, Pole-mounting bracket
2 arms 180°, 750 mm, Black

Recalling the iconic beauty of 20th-century street lanterns, ClassicStreet combines elegant design with the subtle details of heritage luminaires while meeting today's urban lighting requirements. Though echoing the past, the ClassicStreet luminaire is very much shaped for the future. Its gently curved design, together with dedicated brackets and poles, creates a contemporary look that beautifies our cities and creates distinctive and warm night- and daytime ambiances. Not only in historical and traditional architectural settings, but also in more modern urban environments too. Designed around its LED engine, ClassicStreet delivers outstanding levels of lighting and energy performance.

Product data

General Information		Full product name	
Bracket type	Pole-mounting bracket 2 arms 180°	Full product code	JGB795 MBP-T 750 76 BK
Bracket color	Black	Order code	871794321654300
Product family code	JGB795 [CLASSICSTREET BRACKETS POST TOP]	Material Nr. (12NC)	21654300
Value ladder	Specification	Numerator - Quantity Per Pack	912300023377
CE mark	Yes	EAN/UPC - Product/Case	1
Mechanical and Housing		Numerator - Packs per outer box	8717943216543
Bracket length	750 mm	EAN/UPC - Case	1
Product Data			
Order product name	JGB795 MBP-T 750 76 BK		

ClassicStreet brackets

Dimensional drawing

