



DEB Emergency Ballast

Fluorescent Emergency Ballast, 1100-1400 Lumens, One or Two Lamp, T8/T9/T10/T12/Long Compact, Self-Test/Self-Diagnostics

Fluorescent Emergency Ballast, 1100-1400 Lumens, One or Two Lamp, T8/T9/T10/T12/Long Compact, Self-Test/Self-Diagnostics

The DEB Fluorescent emergency ballast carries the load of the luminaire in the event of a power loss, until the generator comes online. They light a specified number of lamps about half brightness to assist in egress from the building.

Product data

Operating and Electrical		Material Nr. (12NC)	912401025072
Voltages	120/277 V	Net Weight (Piece)	1.800 kg
Product Data		Catalog Number	DEB-5ST
Order product name	DEB-5ST	Catalog Number Description	Fluorescent Emergency Ballast, 1100-1400 Lumens, One or Two Lamp, T8/T9/T10/T12/Long Compact, Self-Test/Self-Diagnostics
EAN/UPC - Product	734340860184		
Order code	DEB-5ST		
Numerator - Quantity Per Pack	1		
Numerator - Packs per outer box	1		