



Chronosense

LLC7120/00 CHRONOSENSE

Luminaire based, stand-alone controller for magnetic ballasts up to 600 W. The design of Chronosense is optimized for mounting in a luminaire or gear tray. Step-dim schedule is created using 5 dip switches, dim schedule is flexible between 4 and 9 hours around the midpoint of the night.

Product data

Product Data	
Order product name	LLC7120/00 CHRONOSENSE
Full product name	LLC7120/00 CHRONOSENSE
Full product code	872790086814200
Order code	913700334903
Material Nr. (12NC)	913700334903

Numerator - Quantity Per Pack	1
EAN/UPC - Product/Case	8711500998477
Numerator - Packs per outer box	48
EAN/UPC - Case	8727900868142

Dimensional drawing

