



# Luminaire Based Accessories

## LCA8001/00 RING FOR COVER SET OF 100PCE

Luminaire Based Accessories cover a generic range of products which can be used to extend functionality and/or simplify the design-in of our luminaire based systems.

### Product data

Product Data	
Order product name	LCA8001/00 RING FOR COVER SET OF 100PCE
Full product name	LCA8001/00 RING FOR COVER SET OF 100PCE
Full product code	872790088278000
Local code description	LCA8001/00
Order code	497990
Material Nr. (12NC)	913700338303

Local order code	497990
Numerator - Quantity Per Pack	1
EAN/UPC - Product/Case	8727900882780
Numerator - Packs per outer box	40
EAN/UPC - Case	8727900882797