



Residential Collection

# Glad

## Glad pendant white 1x20W 230V


Glad - pendant - white - 1 pc - 20 W

### Product data

General Information	
Number of light sources	1 pc
Lamp holder	E27
Protection class IEC	Safety class I
Family name	Glad
Product type	pendant
Lamp type	PLCE
Operating and Electrical	
Maximum lamp power	20 W
Approval and Application	
Ingress protection code	IP20 [ Finger-protected]

Product Data	
Full product code	403993116
Order product name	Glad pendant white 1x20W 230V
EAN/UPC - Product	8717943003624
Order code	915002091901
Numerator - Quantity Per Pack	1
Numerator - Packs per outer box	2
Material Nr. (12NC)	915002091901
Net Weight (Piece)	1.998 kg

Dimensional drawing



Residential Collection pendant

