



Residential Collection

Flowerbed

Flowerbed pedestal anthracite 1x11W 230V

Flowerbed - pedestal - antracit - 1 pc - 11 W


Product data

General Information	
Number of light sources	1 pc
Lamp holder	2G7
Protection class IEC	Safety class II
Family name	Flowerbed
Product type	pedestal
Lamp type	PLS4
Operating and Electrical	
Maximum lamp power	11 W
Approval and Application	
Ingress protection code	IP44 [Wire-protected, splash-proof]

Product Data	
Full product code	169159316
Order product name	Flowerbed pedestal anthracite 1x11W 230V
EAN/UPC - Product	8717943747122
Order code	915002543501
Numerator - Quantity Per Pack	1
Numerator - Packs per outer box	3
Material Nr. (12NC)	915002543501
Net Weight (Piece)	1.400 kg

Flowerbed

Dimensional drawing



Residential Collection pedestal

