



Residential Collection

Alive

Alive pendant white 1x20W 230V


Alive - pendant - white - 1 pc - 20 W

Product data

General Information	
Number of light sources	1 pc
Lamp holder	E27
Protection class IEC	Safety class I
Family name	Alive
Product type	pendant
Lamp type	PLCE
Operating and Electrical	
Maximum lamp power	20 W
Approval and Application	
Ingress protection code	IP20 [Finger-protected]

Product Data	
Full product code	408555616
Order product name	Alive pendant white 1x20W 230V
EAN/UPC - Product	8718291462361
Order code	915004217901
Numerator - Quantity Per Pack	1
Numerator - Packs per outer box	2
Material Nr. (12NC)	915004217901
Net Weight (Piece)	1.890 kg

Dimensional drawing



Residential Collection pendant

