



GreenAccent

ST181T A LED20/830 PSU-E MB 24G 1C WH

LED Module 2000 lm, 830 warm white, Power supply unit external (On/Off), Medium beam, White

The new GreenAccent LED accent luminaire has a refined finish, and the 4 different mounting options make it possible to address an environment's accent lighting needs while retaining the same visual concept. With optimized service life in the segment and high energy efficiency, it allows the replacement of traditional metal vapor lamp systems with high quality lighting and faithful color rendering.

Product data

General Information	
Lamp family code	LED20 [LED Module 2000 lm]
Remarks	*-Per Lighting Europe guidance paper "Evaluating performance of LED based luminaires - January 2018": statistically there is no relevant difference in lumen maintenance between B50 and for example B10. Therefore, the median useful life (B50) value also represents the B10 value.
Light Technical	
Number of light sources	1
Light source color	830 warm white
Optic type	Medium beam
Controls and Dimming	
Driver/power unit/transformer	Power supply unit external (On/Off)

Mechanical and Housing	
Housing Color	White
Application Conditions	
Performance ambient temperature	25 °C
Tq	
Product Data	
Order product name	ST181T A LED20/830 PSU-E MB 24G 1C WH
Order code	919406000801
Numerator - Quantity Per Pack	1
Numerator - Packs per outer box	1
Material Nr. (12NC)	919406000801
Full product name	ST181T A LED20/830 PSU-E MB 24G 1C WH