



OmniSpot LED Track Heads

OmniSpot LED 1000lm 90CRI 30K White ELV CrispWhite

OmniSpot LED 1000lm 90CRI 30K White ELV CrispWhite

OmniSpot LED cylinder is designed to meet the most demanding retail challenges through sleek and discrete aesthetics while providing a high center beam candlepower and efficacy performance. OmniSpot LED cylinder allows for less than 1" of space between the track head and track in the ceiling.

Product data

General Information		Initial Performance (IEC Compliant)	
Lamp Source Code	LED	Luminaire Lumens	1000
Colour Temperature	3000		
CRI	90		
Mounting	Track		
Operating and Electrical		Product Data	
Voltages	120 V	Order product name	OmniSpot LED 1000lm 90CRI 30K White ELV CrispWhite
Luminaire Input Watts	12	Order code	LC10930WTECW
		Material Nr. (12NC)	LIGU60014637
		Catalog Number	LC10930WTECW
		Catalog Number Description	OmniSpot LED 1000lm 90CRI 30K White ELV CrispWhite
Mechanical and Housing			
Shape	Round		