



# Calculite LED 4" cylinder gen 3

## Calculite Stem Kit 36" Black

### Calculite Stem Kit 36" Black

Calculite 4" Cylinder Gen3 features industry leading visual comfort and excellent uniform illumination over time.

#### Product data

<b>General Information</b>		<b>Commercial Code</b>	CASK36BK
Lamp Source Code	LED	Order code	CASK36BK
Downlight Type	Cylinders	Catalog Number Description	Calculite Stem Kit 36" Black
Mounting	Suspended	Service / Spec Smart	-
<b>Product Data</b>		Material Nr. (12NC)	LIGU60014711
Order product name	Calculite Stem Kit 36" Black		