

PHILIPS

Ceiling light

myLiving

31141

white

LED



31141/31/66



Light and relaxed


Accentuate your home with the versatile Philips myLiving Tigo LED wall and ceiling light White. Its sand blasted glass cover secured with robust metal clips, provides diffused warm white light for that cosy feel.

Pure, natural light

- Warm white light effect

Highlights

Warm white light effect



Light up your home with warm white light effects thanks to the latest LED innovation. Philips introduces the latest generation of warm white light high power LEDs, reaching a colour temperature of 2700 Kelvin. This innovation brings the Philips lighting fixtures to the next level, allowing you to combine elegance and ambience

creation with light and energy efficiency.

Specifications

Design and finishing

- Color: white
- Material: glass

Miscellaneous

- Especially designed for: Living- & Bedroom
- Style: Contemporary
- Type: Ceiling light
- EyeComfort: No

Product dimensions & weight

- Height: 40 cm
- Length: 40 cm
- Net weight: 1.420 kg
- Width: 12.7 cm

Service

- Warranty: 1 year(s)

Technical specifications

- Mains power: 230 V
- Fixture dimmable: No
- LED: Yes
- IP code: IP20, protection against objects bigger than 12.5 mm, no protection against water
- Class of protection: I - earthed



Issue date: 2023-11-08
Version: 0.487

© 2023 Signify Holding. All rights reserved. Signify does not give any representation or warranty as to the accuracy or completeness of the information included herein and shall not be liable for any action in reliance thereon. The information presented in this document is not intended as any commercial offer and does not form part of any quotation or contract, unless otherwise agreed by Signify. Philips and the Philips Shield Emblem are registered trademarks of Koninklijke Philips N.V.

www.lighting.philips.com