



Renaissance Bed Light

4' Economy Patient Bed Light, 2-lamps up/2-lamps down, 9618 Series

4' Economy Patient Bed Light, 2-lamps up/2-lamps down, 9618 Series

Economy patient bed light with separately switched direct/indirect illumination. Prismatic acrylic diffusers provide quality illumination.

Product data

General Information	
Lamp Source Code	FLU
Lamp Watts	32
Lamp Envelope Code	T-8
Number of Lamps	4
Photometric Distribution Type	General Diffuse
Operating and Electrical	
Voltages	120/277/347 V
Mechanical and Housing	
Shape	Linear

Product Data	
Order product name	4' Economy Patient Bed Light, 2-lamps up/2-lamps down, 9618 Series
Commercial Code	9618-422
Order code	9618-422
Catalog Number Description	4' Economy Patient Bed Light, 2-lamps up/2-lamps down, 9618 Series
Material Nr. (12NC)	ALK00000849