



Multi-Med Recessed

Multi-Med 2x4 Recessed Exam/Ambient Patient Room Luminaire

Multi-Med 2x4 Recessed Exam/Ambient Patient Room Luminaire

The Multi-Med Patient Luminaire provides the three essential illumination requirements at the patient bed – ambient, examination and optional night/chart lighting. The optical compartments are designed for efficient performance and patient comfort. The overall clean appearance of Multi-Med allows it to blend into the ceiling plane – helping the environment to feel less clinical. The standard antimicrobial finish and optional fully-lensed door model supports a healthcare facility's infection-control efforts.

Product data

General Information		Shape	
Lamp Watts	40	Rectangular	
Lamp Envelope Code	TT5		
Number of Lamps	5		
Photometric Distribution Type	Direct		
Operating and Electrical		Product Data	
Voltages	120/277/347 V	Order product name	Multi-Med 2x4 Recessed Exam/Ambient Patient Room Luminaire
		Commercial Code	ME540N-120-PBL-(AMBIENT)
		Order code	ME540N-120-PBL-(AMBIENT)
		Catalog Number Description	Multi-Med 2x4 Recessed Exam/Ambient Patient Room Luminaire
		Service / Spec Smart	-
		Material Nr. (12NC)	ALK00004032
Mechanical and Housing			
Fixture Size	2 ft x 4 ft		