



# Multi-Med Recessed

## Multi-Med 2x2 Recessed Exam Patient Room Luminaire

### Multi-Med 2x2 Recessed Exam Patient Room Luminaire

The Multi-Med Patient Luminaire provides the three essential illumination requirements at the patient bed – ambient, examination and optional night/chart lighting. The optical compartments are designed for efficient performance and patient comfort. The overall clean appearance of Multi-Med allows it to blend into the ceiling plane – helping the environment to feel less clinical. The standard antimicrobial finish and optional fully-lensed door model supports a healthcare facility's infection-control efforts.

#### Product data

General Information	
Lamp Watts	40
Lamp Envelope Code	TT5
Number of Lamps	4
Operating and Electrical	
Voltages	120/277/347 V
Mechanical and Housing	
Fixture Size	2 ft x 2 ft
Shape	Square

Product Data	
Order product name	Multi-Med 2x2 Recessed Exam Patient Room Luminaire
Commercial Code	ME440NLG
Order code	ME440NLG
Catalog Number Description	Multi-Med 2x2 Recessed Exam Patient Room Luminaire
Service / Spec Smart	-
Material Nr. (12NC)	ALK00004042