



SurgiCare Recessed 2x4

2x4 Supplemental Surgical Luminaire, STR24 Series

2x4 Supplemental Surgical Luminaire, STR24 Series

The STR Series is designed for the illumination, electrical shielding and infection control requirements of the surgical suite. The asymmetric light distribution provides the ambient illumination while contributing supplemental illumination for the surgical field. STR luminaires are available in a variety of lamping options. Meets MIL-STD 461F.

Product data

General Information	
Lamp Watts	28
Lamp Envelope Code	T-5
Number of Lamps	6
Operating and Electrical	
Voltages	120/277 V
Mechanical and Housing	
Fixture Size	2 ft x 4 ft
Shape	Rectangular

Product Data	
Order product name	2x4 Supplemental Surgical Luminaire, STR24 Series
Commercial Code	STR(F/G)24-6-T5
Order code	STR(F/G)24-6-T5
Catalog Number Description	2x4 Supplemental Surgical Luminaire, STR24 Series
Service / Spec Smart	-
Material Nr. (12NC)	ALK00004139