



SurgiCare Recessed Sealed 2x4

2x4 Sealed Supplemental Surgical Luminaire, SSTR24 Series

2x4 Sealed Supplemental Surgical Luminaire, SSTR24 Series

The SSTR Series is designed for the illumination, electrical shielding and infection control requirements of the surgical suite. The asymmetric/symmetric light distribution provides the ambient illumination while contributing supplemental illumination for the surgical field.

SSTR luminaires are available in a variety of lamping options. The housing is fully sealed to minimize potential room contamination and simplify cleaning. Meets MIL-STD 461F, wet location listed and is IP65 rated.

Product data

General Information	
Lamp Source Code	FLU
Lamp Watts	28
Lamp Envelope Code	T-5
Number of Lamps	3
Operating and Electrical	
Voltages	120/277 V
Mechanical and Housing	
Fixture Size	2 ft x 4 ft

Shape	Rectangular
Product Data	
Order product name	2x4 Sealed Supplemental Surgical Luminaire, SSTR24 Series
Commercial Code	SSTR(F/G)24-3-T5
Order code	SSTR(F/G)24-3-T5
Catalog Number Description	2x4 Sealed Supplemental Surgical Luminaire, SSTR24 Series
Service / Spec Smart	-
Material Nr. (12NC)	ALK00004170