



DormLyte

Dorm Lyte Wall Mounted Fluorescent, 6" Width, 48" Nominal Length containing one 32 watt T8 lamp

Dorm Lyte Wall Mounted Fluorescent, 6" Width, 48" Nominal Length
- containing one 32 watt T8 lamp

Dormlyte Wall Mount: For dorms, nursing homes, assisted living, health care facilities, or other over bed or wall mounted applications. 1/8" thick DR acrylic lens. Unique co-extension. White opaque center. Frosted white front sections for subtle glow. Linear prisms to spread up/downlight. Only 3-7/8" deep, meets American With Disabilities Act for wall mounted fixtures. Multiple lumen/watt packages for design flexibility and added energy savings.

Product data

General Information	
Lamp Source Code	FLU
Lamp Watts	32
Lamp Envelope Code	T-8
Number of Lamps	1
Mounting	Wall Surface
Operating and Electrical	
Voltages	277 V
Mechanical and Housing	
Fixture Size	48 in

Product Data	
Order product name	Dorm Lyte Wall Mounted Fluorescent, 6" Width, 48" Nominal Length containing one 32 watt T8 lamp
Commercial Code	DM4X132UNVP2
Order code	DM4X132UNVP2
Catalog Number Description	Dorm Lyte Wall Mounted Fluorescent, 6" Width, 48" Nominal Length containing one 32 watt T8 lamp
Material Nr. (12NC)	CFI00001499