



Product Description

MASTER PL-T 2 Pin

Energy saving compact fluorescent lamp Compact long-arc low-pressure mercury discharge lamp Envelope consists of 6 parallel narrow fluorescent tubes

Benefits

- Optimal light performance for short or long operation per switch on
- Substantial more light (about 15% more) per light point than common
- Long life and sustainable light output

Features

- 2-pin base with a housing incorporating the lamp glow switch starter
- Fast light run up
- Philips invented bridge technology applied
- Dimming is not possible

Application


- Designed for general or supplementary lighting in professional and semi-professional environments
- Ideal for use in recessed or surface mounted downlighters, wall units

Versions



MASTER PL-T 2 Pin

Dimensional drawing



Product	D (max)	A (max)	B (max)	C (max)
MASTER PL-T 26W/840/2P 1CT/ 5X10BOX	41 mm	86.7 mm	111.0 mm	134.2 mm
MASTER PL-T 26W/830/2P 1CT/ 5X10BOX	41 mm	86.7 mm	111.0 mm	134.2 mm

Approval and Application

Energy Consumption kWh/1000 h	32 kWh
Energy Efficiency Label (EEL)	B
Mercury (Hg) Content (Nom)	1.4 mg

Controls and Dimming

Dimmable	no
----------	----

Operating and Electrical

Lamp Current (Nom)	0.335 A
Power (Rated) (Nom)	26.0 W

General Information

Cap-Base	GX24D-3
Life To 10% Failures (Nom)	6500 h
Life To 50% Failures (Nom)	10000 h
LSF 2000 h Rated	99 %
LSF 4000 h Rated	98 %
LSF 6000 h Rated	92 %
LSF 8000 h Rated	78 %

Light Technical

LLMF 2000 h Rated	92 %
LLMF 4000 h Rated	87 %
LLMF 6000 h Rated	84 %
LLMF 8000 h Rated	81 %
Luminous Efficacy (rated) (Min)	66 lm/W
Luminous Flux (Nom)	1800 lm
Luminous Flux (Rated) (Nom)	1725 lm

Temperature

Design Temperature (Nom)	28 °C
--------------------------	-------

Light Technical

Order Code	Full Product Name	Colour Code	Colour Designation	Correlated Colour Temperature (Nom)	Colour Rendering
					Index (Nom)
927914583071	MASTER PL-T 26W/ 830/2P 1CT/5X10BOX	830	Warm white (WW)	3000 K	>80
927914584071	MASTER PL-T 26W/ 840/2P 1CT/5X10BOX	840	Cool White (CW)	4000 K	82

MASTER PL-T 2 Pin

