



Derungs Bed Light – stylish multi-purpose solution

Derungs, room light

The Derungs Bed Light is an elegant, versatile solution offering a choice of four settings: night light, general light (up), examination light (all on) and reading light. This one luminaire can therefore cover all the functions required in a patient room. Glare-free prisms ensure a comfortable lighting effect.

Benefits

- Multi-functional – one luminaire performs all functions needed in a patient room
- Warm, gentle glare-free light
- Durable, stylish design

Features

- Four light sources providing all required settings (can be switched separately)
- Compliant with relevant norms, including EN 60598
- DRS double-reflecting system for glare-free and efficient lighting
- Night-light included

Application

- Patient rooms

Specifications

Type	CWS490
------	--------

Light source	Fluorescent:
	- 2 x MASTER TL5 / G5 / 39 W
	- 1 x MASTER TL5 / G5 / 24 W
	- 1x GLS / E14 / MAX 7 W

Derungs, room light

Lamp included	Yes (K)
Ballast	Electronic ballast
Connection	Mains lead with special plug, approx 2 m
Material	Luminaire body: aluminum, powder coated
	Lamp cover: acrylic satine


Installation	Screw mounting, directly onto wall
	Operating temperature: 0 °C to +40 °C
Remarks	EN 60598-2-25, EN 60598-1

Versions



Derungs CWS490 room light

Dimensional drawing



DuraCare CWS/TWS490/491

