



Aesthetically pleasing and great flexibility

DuoCare Recessed Exam/Ambient

DuoCare Exam/Ambient Patient Room Luminaire provides ambient and examination illumination for the patient room. The two luminaire system utilizes an asymmetric light distribution to reduce glare for the patient and shadows for the healthcare professional. A narrow aperture and low housing height provides greater flexibility in ceiling placement.

Benefits

- Easy relamping

Features

- Housing and door frame have an antimicrobial paint finish

Application

- Hospitals and Care Centers

Versions

