



Helps the environment feel less medical

Multi-Med Recessed

The Multi-Med Patient Luminaire provides the three essential illumination requirements at the patient bed – ambient, examination and optional night/chart lighting. The optical compartments are designed for efficient performance and patient comfort. The overall clean appearance of Multi-Med allows it to blend into the ceiling plane – helping the environment to feel less clinical. The standard antimicrobial finish and optional fully-lensed door model supports a healthcare facility's infection-control efforts.

Benefits

- Tool-less access for lamp replacement

Features

- Antimicrobial finish

Application

- Hospitals and Care Centers

Multi-Med Recessed

Versions

