



# Meets all illumination needs for a surgical suite

## SurgiCare Recessed 2x4

The STR Series is designed for the illumination, electrical shielding and infection control requirements of the surgical suite. The asymmetric light distribution provides the ambient illumination while contributing supplemental illumination for the surgical field. STR luminaires are available in a variety of lamping options. Meets MIL-STD 461F.

### Benefits

- One piece door frame
- antimicrobial paint finish
- smooth, room side lens

### Features

- Meets MIL-STD 461F

### Application

- Hospitals and Care Centers

### Versions

