



eW Profile Powercore

ZCX411 RA L483 (25PCS)

Mounting Track - 483 mm (25 @ 19 in) - White

Mounting Track - 483 mm (25 @ 19 in) - White

Product data

General Information	
Type description	MOUNTING KIT [Mounting kit]
Suspension accessories	RA [Mounting rail/profile]
Packaging configuration	(25PCS) [25 pieces]
Accessory color	White

Operating and Electrical	
Electrical accessories	No [-]
Input Voltage	230 V

Mechanical and Housing	
Length	483 mm

Product Data	
Full product code	872790071474600
Order product name	ZCX411 RA L483 (25PCS)
EAN/UPC - Product	8727900714746
Order code	910503700271
Numerator - Quantity Per Pack	1
Numerator - Packs per outer box	1
Material Nr. (12NC)	910503700271
Net Weight (Piece)	3.570 kg