



StyliD PerfectBeam

ST741Z LBSP

StyliD PerfectBeam has been developed specifically for museum applications. It strikes the perfect balance between presenting art and preserving it, between ever-changing exhibitions and the need for low maintenance and energy efficiency. The often conflicting demands for good visibility and conservation are now being met with LED lighting that emits no harmful UV or IR radiation. StyliD PerfectBeam presents works of art in their true colors, subtly attracting the viewer's eye, making it easy to see, interpret and appreciate masterpieces down to the smallest detail. It enables you to use light to evoke emotions and create an unforgettable museum experience. Its modular design, variable beam, beam shaping accessories and onboard dimming ensure maximum flexibility to cope with changing exhibitions and lighting needs, now and in the future.

Product data

| Mechanical and Housing | | Order code | |
|-------------------------------------|---|---------------------------------|--------------|
| Dimensions (height x width x depth) | NaN x NaN x NaN mm (NaN x NaN x NaN in) | 29783300 | |
| Product Data | | Numerator – quantity per pack | 1 |
| Full product code | 871829129783300 | Numerator – packs per outer box | 1 |
| Order product name | ST741Z LBSP | Material no. (12NC) | 910504801189 |
| EAN/UPC – product | 8718291297833 | Net weight (piece) | 18.000 g |