



Latina FBH024

GBH024 FRG

Latina is a recessed downlight for compact fluorescent PL-C lamps, with spring fasteners for easy and quick installation. Square and round models are available. Latina comes in ready-to-install kit versions, including lamps and pre-connected separate gearbox. The round models are available in glass versions (frosted glass or clear decorative glass) or in open versions. For the open versions, the same glass designs are available as accessories (to be ordered separately). The square versions are always with glass.

Product data

Product Data	
Order product name	GBH024 FRG
Full product name	GBH024 FRG
Full product code	871794388531299
Order code	910930025218
Material Nr. (12NC)	910930025218

Numerator - Quantity Per Pack	1
EAN/UPC - Product/Case	8717943885312
Numerator - Packs per outer box	20
EAN/UPC - Case	8717943886753