



ClassicStreet brackets

JGB794 MBP-T 750 76 BK

CLASSICSTREET BRACKETS SUSPENDED, Pole-mounting bracket 2 arms 180°, 750 mm, Black

Recalling the iconic beauty of 20th-century street lanterns, ClassicStreet combines elegant design with the subtle details of heritage luminaires while meeting today's urban lighting requirements. Though echoing the past, the ClassicStreet luminaire is very much shaped for the future. Its gently curved design, together with dedicated brackets and poles, creates a contemporary look that beautifies our cities and creates distinctive and warm night- and daytime ambiances. Not only in historical and traditional architectural settings, but also in more modern urban environments too. Designed around its LED engine, ClassicStreet delivers outstanding levels of lighting and energy performance.

Product data

General Information			
Bracket type	Pole-mounting bracket 2 arms 180°	Full product name	JGB794 MBP-T 750 76 BK
Bracket color	Black	Full product code	871794321652900
Product family code	JGB794 [CLASSICSTREET BRACKETS SUSPENDED]	Order code	912300023375
Value ladder	Specification	Material Nr. (12NC)	912300023375
CE mark	Yes	Numerator - Quantity Per Pack	1
Mechanical and Housing		EAN/UPC - Product/Case	8717943216529
Bracket length	750 mm	Numerator - Packs per outer box	1
Product Data		EAN/UPC - Case	8717943216529
Order product name	JGB794 MBP-T 750 76 BK		

ClassicStreet brackets

Dimensional drawing

