



TownTune Central Post-Top DR

ZDP261 DR GO

TOWNTUNE ACC., 5 J vandal-protected, Gold

Designed to enhance existing and scalable urban spaces, the Philips TownTune family offers all the latest lighting innovations in terms of performance, quality of light and connectivity. The family consists of four solutions: a Central Post Top (CPT), an Asymmetric Spigot Post Top / Side Entry version (ASY), a version with an extending Lyre post top bracket (Lyre), and a Central Post Top with a Conical Comfort Bowl (CCB). Each TownTune luminaire can be customised with a choice of different shapes on top of the housing, plus there's the option to add a decorative ring, which comes in two colours (excluding CCB). Design options that enable you to create your very own lighting signature and bring a distinctive identity to districts and cities. In addition, every luminaire in the TownTune family is uniquely identifiable, thanks to the Signify Service tag app. By simply scanning a QR code, placed inside the door of the mast or directly on the luminaire, you can instantly access the configuration of the luminaire. This makes maintenance and programming operations faster and easier and enables you to create your digital library of lighting assets and spare parts. TownTune also uses the Philips LEDGINE-O lighting platform, ensuring you always have the right amount and direction of light on your street. Furthermore, thanks to being system ready (SR), TownTune is also future proof. A solution that's ready to be paired with both standalone and advanced control and lighting software applications such as Interact City.

Product data

General Information		Product family code	ZDP261 [TOWNTUNE ACC.]
Angle	-	Value ladder	Performance
Accessory colour	Gold		

TownTune Central Post-Top DR

Mechanical and Housing	
Overall height	61 mm
Overall diameter	533 mm
Material	Polymethyl methacrylate
Approval and Application	
Ingress protection code	-
Mech. impact protection code	IK08 [5 J vandal-protected]
Product Data	
Order product name	ZDP261 DR GO

Full product name	ZDP261 DR GO
Full EOC	871869949033100
Order code	49033100
Material no. (12 NC)	912300024164
SAP numerator – quantity per pack	1
EAN/UPC – Product/Case	8718699490331
Numerator – packs per outer box	1
EAN/UPC – Case	8718699490331

