



Blast gen4

ZCP485 TRIM RING WH

Diffuser Trim Ring - White


Diffuser Trim Ring - White

Product data

General Information	
Accessory color	White
Product Data	
Full product code	871829138812899
Order product name	ZCP485 TRIM RING WH
EAN/UPC - Product	8718291388128
Order code	912400130336

Local order code	.
Numerator - Quantity Per Pack	1
Numerator - Packs per outer box	20
Material Nr. (12NC)	912400130336
Net Weight (Piece)	0.290 kg

Dimensional drawing



Blast gen4 acc. ZCP845