



MTD Series Auxillary Circuit Transfer Device

MTD

Auxiliary Circuit Transfer Device, For Internal Luminaire Mounting

The MTD Auxiliary Power Transfer Device is designed to transfer generator or AC systems power to AC fluorescent lighting ballasts in the event of a utility power failure. This transfer function occurs independent of wall switch position.

Product data

Operating and Electrical		Numerator - Packs per outer box 10	
Voltages	120/277 V	Material Nr. (12NC)	912400370843
Product Data		Net Weight (Piece)	0.001 kg
Order product name	MTD	Commercial Code	MTD
EAN/UPC - Product	781370584191	Catalog Number Description	Auxiliary Circuit Transfer Device, For Internal Luminaire Mounting
Order code	912400370843	Quick Ship	Service Smart (2 Days) USA
Local order code	7058419		
Numerator - Quantity Per Pack	1		