



Feroe Brackets

JRP531 MBP-T 1000 60 RAL

Pole-mounting bracket 2 arms 180°, 60 mm, RAL color

The Feroe bracket displays taut lines reinforced by a lining on the inside that accentuate its finesse and command at the same time.

Product data

General Information	
Bracket type	Pole-mounting bracket 2 arms 180°
Bracket color	RAL color
Value ladder	Performance
CE mark	Yes
Mechanical and Housing	
Fixation point length distance (Nom)	100 mm
Fixation diameter	Mounting device diameter 60 mm
Effective projected area	0.094 m²
Bracket tilt angle	5°
Bracket to luminaire diameter	60 mm
Bracket length	1000 mm

Product Data	
Order product name	JRP531 MBP-T 1000 60 RAL
Full product name	JRP531 MBP-T 1000 60 RAL
Full product code	871869693278000
Order code	912401451171
Material Nr. (12NC)	912401451171
Numerator - Quantity Per Pack	1
EAN/UPC - Product/Case	8718696932780
Numerator - Packs per outer box	1
EAN/UPC - Case	8718696932780

Feroe Brackets

Dimensional drawing

