



Feroe Brackets

JRP530 MBP-T 1000 76 RAL

Pole-mounting bracket 2 arms 180°, 76 mm, RAL color

The Feroe bracket displays taut lines reinforced by a lining on the inside that accentuate its finesse and command at the same time.

Product data

General Information	
Bracket type	Pole-mounting bracket 2 arms 180°
Bracket color	RAL color
Value ladder	Performance
CE mark	Yes
Mechanical and Housing	
Fixation point length distance (Nom)	100 mm
Fixation diameter	Mounting device diameter 60 mm
Effective projected area	0.12 m²
Bracket tilt angle	5°
Bracket to luminaire diameter	76 mm
Bracket length	1000 mm

Product Data	
Order product name	JRP530 MBP-T 1000 76 RAL
Full product name	JRP530 MBP-T 1000 76 RAL
Full product code	871869693285800
Order code	912401451178
Material Nr. (12NC)	912401451178
Numerator - Quantity Per Pack	1
EAN/UPC - Product/Case	8718696932858
Numerator - Packs per outer box	1
EAN/UPC - Case	8718696932858

Feroe Brackets

Dimensional drawing

