



# KC Brackets

## JGM001 KC MBP-T 500 60 RAL

Pole-mounting bracket 2 arms 180°, 500 mm, 60 mm, 5°, RAL color, Mounting device diameter 60 mm

The KC bracket is functional, adjustable and will fit any type of landscape with the utmost discretion. Its minimalistic design is timeless. Furthermore, five different bracket lengths allow various configurations to provide efficient lighting through any urban constraint.

### Product data

| General Information                  |                                   |
|--------------------------------------|-----------------------------------|
| Bracket type                         | Pole-mounting bracket 2 arms 180° |
| Bracket color                        | RAL color                         |
| Value ladder                         | Value                             |
| CE mark                              | Yes                               |
| Mechanical and Housing               |                                   |
| Fixation point length distance (Nom) | 70 mm                             |
| Fixation diameter                    | Mounting device diameter 60 mm    |
| Effective projected area             | 0.082 m²                          |
| Bracket tilt angle                   | 5°                                |
| Bracket to luminaire diameter        | 60 mm                             |
| Bracket length                       | 500 mm                            |

| Product Data                    |                            |
|---------------------------------|----------------------------|
| Order product name              | JGM001 KC MBP-T 500 60 RAL |
| Full product name               | JGM001 KC MBP-T 500 60 RAL |
| Full product code               | 871869693310700            |
| Order code                      | 912401451204               |
| Material Nr. (12NC)             | 912401451204               |
| Numerator - Quantity Per Pack   | 1                          |
| EAN/UPC - Product/Case          | 8718696933107              |
| Numerator - Packs per outer box | 1                          |
| EAN/UPC - Case                  | 8718696933107              |

## KC Brackets

### Dimensional drawing

