



Tradition Brackets

JSP299 MBP 1000 60 RAL

Pole-mounting bracket, 1000 mm, 60 mm, 0°, RAL color, Mounting device diameter 1"

The Tradition bracket has a long-established and simple shape which powdered-steel finish gives an up-to-date and relevant presence in today's landscapes.

Product data

General Information	
Bracket type	Pole-mounting bracket
Bracket color	RAL color
Value ladder	Performance
CE mark	Yes
Mechanical and Housing	
Fixation point length distance (Nom)	100 mm
Fixation diameter	Mounting device diameter 1"
Effective projected area	0.23 m²
Bracket tilt angle	0°
Bracket to luminaire diameter	60 mm
Bracket length	1000 mm

Product Data	
Order product name	JSP299 MBP 1000 60 RAL
Full product name	JSP299 MBP 1000 60 RAL
Full product code	871869693374900
Order code	912401451268
Material Nr. (12NC)	912401451268
Numerator - Quantity Per Pack	1
EAN/UPC - Product/Case	8718696933749
Numerator - Packs per outer box	1
EAN/UPC - Case	8718696933749

Tradition Brackets

Dimensional drawing

