



Wy Brackets

JGB700 MBP-T 1000 60 RAL

Pole-mounting bracket 2 arms 180°, 1000 mm, 60 mm, 0°, RAL color, Mounting device diameter 3/4"

The Wy bracket enriches a landscape with elegance and sophistication thanks to some very elaborated details.

Product data

General Information	
Bracket type	Pole-mounting bracket 2 arms 180°
Bracket color	RAL color
Value ladder	Performance
CE mark	Yes
Mechanical and Housing	
Fixation point length distance (Nom)	25 mm
Fixation diameter	Mounting device diameter 3/4"
Effective projected area	0.13 m²
Bracket tilt angle	0°
Bracket to luminaire diameter	60 mm
Bracket length	1000 mm

Product Data	
Order product name	JGB700 MBP-T 1000 60 RAL
Full product name	JGB700 MBP-T 1000 60 RAL
Full product code	871869693394700
Order code	912401451288
Material Nr. (12NC)	912401451288
Numerator - Quantity Per Pack	1
EAN/UPC - Product/Case	8718696933947
Numerator - Packs per outer box	1
EAN/UPC - Case	8718696933947

Wy Brackets

Dimensional drawing

