



Notre Dame Brackets

JGP300 MBP-S 700 89 RAL

Pole-mounting bracket 1 arm, 700 mm, 89 mm, 0°, RAL color, Mounting device diameter 1"

The Notre Dame bracket makes a strong statement in a landscape, mainly highlighting the authenticity of the surroundings and giving them some contemporary edge.

Product data

General Information	
Bracket type	Pole-mounting bracket 1 arm
Bracket color	RAL color
Value ladder	Performance
CE mark	Yes
Mechanical and Housing	
Fixation point length distance (Nom)	45 mm
Fixation diameter	Mounting device diameter 1"
Effective projected area	0.2621 m²
Bracket tilt angle	0°
Bracket to luminaire diameter	89 mm
Bracket length	700 mm

Product Data	
Order product name	JGP300 MBP-S 700 89 RAL
Full product name	JGP300 MBP-S 700 89 RAL
Full product code	871869693668900
Order code	912401451456
Material Nr. (12NC)	912401451456
Numerator - Quantity Per Pack	1
EAN/UPC - Product/Case	8718696936689
Numerator - Packs per outer box	1
EAN/UPC - Case	8718696936689

Notre Dame Brackets

Dimensional drawing

