



Notre Dame Brackets

JGP300 MBP-T 700 89 RAL

Pole-mounting bracket 2 arms 180°, 700 mm, 89 mm, 0°, RAL color, Mounting device diameter 1"

The Notre Dame bracket makes a strong statement in a landscape, mainly highlighting the authenticity of the surroundings and giving them some contemporary edge.

Product data

General Information	
Bracket type	Pole-mounting bracket 2 arms 180°
Bracket color	RAL color
Value ladder	Performance
CE mark	Yes
Mechanical and Housing	
Fixation point length distance (Nom)	45 mm
Fixation diameter	Mounting device diameter 1"
Effective projected area	0.3141 m²
Bracket tilt angle	0°
Bracket to luminaire diameter	89 mm
Bracket length	700 mm

Product Data	
Full product code	871869693669600
Order product name	JGP300 MBP-T 700 89 RAL
Order code	912401451457
Numerator - Quantity Per Pack	1
Numerator - Packs per outer box	1
Material Nr. (12NC)	912401451457
Full product name	JGP300 MBP-T 700 89 RAL
EAN/UPC - Case	8718696936696
EAN/UPC - Product/Case	8718696936696

Notre Dame Brackets

Dimensional drawing

