



Notre Dame Brackets

JGP300 MBW 700 RAL

Wall-mounting bracket, 700 mm, 89 mm, 0°, RAL color,
Mounting device diameter 1"

The Notre Dame bracket makes a strong statement in a landscape, mainly highlighting the authenticity of the surroundings and giving them some contemporary edge.

Product data

General Information	
Bracket type	Wall-mounting bracket
Bracket color	RAL color
Value ladder	Performance
CE mark	Yes
Mechanical and Housing	
Fixation point length distance (Nom)	45 mm
Fixation diameter	Mounting device diameter 1"
Effective projected area	0 m²
Bracket tilt angle	0°
Bracket to luminaire diameter	89 mm
Bracket length	700 mm

Product Data	
Full product code	871869693670200
Order product name	JGP300 MBW 700 RAL
Order code	912401451458
Numerator - Quantity Per Pack	1
Numerator - Packs per outer box	1
Material Nr. (12NC)	912401451458
Full product name	JGP300 MBW 700 RAL
EAN/UPC - Case	8718696936702
EAN/UPC - Product/Case	8718696936702

Notre Dame Brackets

Dimensional drawing

