



Notre Dame Brackets

JGP300 MBP-S 1100 89 RAL

Pole-mounting bracket 1 arm, 1100 mm, 89 mm, 0°, RAL color, Mounting device diameter 1"

The Notre Dame bracket makes a strong statement in a landscape, mainly highlighting the authenticity of the surroundings and giving them some contemporary edge.

Product data

General Information	
Bracket type	Pole-mounting bracket 1 arm
Bracket color	RAL color
Value ladder	Performance
CE mark	Yes
Mechanical and Housing	
Fixation point length distance (Nom)	87 mm
Fixation diameter	Mounting device diameter 1"
Effective projected area	0.3341 m²
Bracket tilt angle	0°
Bracket to luminaire diameter	89 mm
Bracket length	1100 mm

Product Data	
Order product name	JGP300 MBP-S 1100 89 RAL
Full product name	JGP300 MBP-S 1100 89 RAL
Full product code	871869693671900
Order code	912401451459
Material Nr. (12NC)	912401451459
Numerator - Quantity Per Pack	1
EAN/UPC - Product/Case	8718696936719
Numerator - Packs per outer box	1
EAN/UPC - Case	8718696936719

Notre Dame Brackets

Dimensional drawing

