



# Notre Dame Brackets

## JGP300 MBP-T 1100 89 RAL

Pole-mounting bracket 2 arms 180°, 1100 mm, 89 mm, 0°, RAL color, Mounting device diameter 1"

The Notre Dame bracket makes a strong statement in a landscape, mainly highlighting the authenticity of the surroundings and giving them some contemporary edge.

### Product data

General Information	
Bracket type	Pole-mounting bracket 2 arms 180°
Bracket color	RAL color
Value ladder	Performance
CE mark	Yes
Mechanical and Housing	
Fixation point length distance (Nom)	87 mm
Fixation diameter	Mounting device diameter 1"
Effective projected area	0.4891 m²
Bracket tilt angle	0°
Bracket to luminaire diameter	89 mm
Bracket length	1100 mm

Product Data	
Order product name	JGP300 MBP-T 1100 89 RAL
Full product name	JGP300 MBP-T 1100 89 RAL
Full product code	871869693672600
Order code	912401451460
Material Nr. (12NC)	912401451460
Numerator - Quantity Per Pack	1
EAN/UPC - Product/Case	8718696936726
Numerator - Packs per outer box	1
EAN/UPC - Case	8718696936726

## Notre Dame Brackets

### Dimensional drawing

