



# OccuSwitch DALI

## LCC2070/00 Cable Basic

Wieland cable for OccuSwitch DALI LRM2070

Wieland cable for OccuSwitch DALI LRM2070

### Product data

Product Data	
Full product code	871155973248099
Order product name	LCC2070/00 CABLE BASIC
EAN/UPC – product	8711559732480
Order code	73248099
Numerator – quantity per pack	1

Numerator – packs per outer box	18
Material no. (12NC)	913700333703
Net weight (piece)	0.100 kg