



OccuSwitch DALI

LCC2080/00 CABLE ADV + BMS

Wieland cable for OccuSwitch DALI LRM2080 and LRM2090

Wieland cable for OccuSwitch DALI LRM2080 and LRM2090

Product data

Product Data	
Full product code	871155973250399
Order product name	LCC2080/00 CABLE ADV + BMS
EAN/UPC - Product	8711559732503
Order code	913700333803
Numerator - Quantity Per Pack	1

Numerator - Packs per outer box	18
Material Nr. (12NC)	913700333803
Net Weight (Piece)	0.100 kg